

Quality Control Certificate



Phenylalanine (PKU) neonatal Screening Assay

REF RE80015 & RE80019	LOT NP258	2026-06-30 YYYY-MM-DD
------------------------------	------------------	---------------------------------

Final Release Results:

	Concentration mg/dL	Measured value OD	Range OD	OD/ODmax
CAL A	1.80	0.052	-	17%
CAL B	3.30	0.084	-	28%
CAL C	6.50	0.118	-	39%
CAL D	11.50	0.216	-	72%
CAL E	17.00	0.300	-	100%

	Concentration found mg/dL	Acceptable range		
		lower limit mg/dL	target value mg/dL	upper limit mg/dL
CONTROL 1	3.1	2.4	3.1	3.7
CONTROL 2	5.1	4.2	5.2	6.2

Kit Composition:

REF		LOT	
KENP655	CONTROL 1	289	2026-06
KENP665	CONTROL 2	289	2026-06
KENP545	LYO ENZ	287	2026-12
KENP755/9	ENZDIL	274	2026-09
KENP615	CAL A-E	288	2026-06
KENP775/9	DYE	275	2026-09
KENP765/9	SUBS	276	2026-09

Valid edition of the instructions for use: 2024-11

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2025-12-16

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl