

# Quality Control Certificate



## alpha-Amylase Saliva Assay

REF RE80111

LOT EAA135

2026-06-30  
YYYY-MM-DD

### Final Release Results:

	Concentration U/mL	Measured value		Delta OD	Acceptabel range
		OD <sub>3min</sub>	OD <sub>8min</sub>		OD <sub>8min</sub>
CAL A	380	1.131	2.675	1.545	< 2.900
CAL B	190	0.582	1.376	0.795	-
CAL C	95	0.310	0.721	0.412	-
CAL D	47.5	0.178	0.399	0.222	-
CAL E	0	0.026	0.026	0.000	-

	Concentration found U/mL	Acceptable range		
		lower limit U/mL	target value U/mL	upper limit U/mL
CONTROL 1	42.3	27.4	45.7	63.9
CONTROL 2	127	86	135	183

### Kit Composition:

REF		LOT	
KEAA651	LYO CONTROL 1	L467	2026-06
KEAA661	LYO CONTROL 2	L468	2026-06
KEAA511	MTP	L469	2026-08
KEAA611	LYO CAL	L466	2026-06
KEAA751	CONC DILBUF	L470	2026-08
KEAA771	SUBS	L471	2026-06

Valid edition of the instructions for use: 2014-04

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

# Accepted

Validation date: 2024-08-06

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)