

Quality Control Certificate



alpha-Amylase Saliva Assay

REF RE80111

LOT EAA136

2027-08-31
YYYY-MM-DD

Final Release Results:

	Concentration U/mL	Measured value		Delta OD	Acceptabel range
		OD _{3min}	OD _{8min}		OD _{8min}
CAL A	340	0.947	2.290	1.343	< 2.900
CAL B	170	0.491	1.197	0.706	-
CAL C	85	0.261	0.627	0.366	-
CAL D	42.5	0.143	0.329	0.185	-
CAL E	0	0.085	0.180	0.095	-

	Concentration found U/mL	Acceptable range		
		lower limit U/mL	target value U/mL	upper limit U/mL
CONTROL 1	41.1	24.1	40.2	56.3
CONTROL 2	135	85	132	180

Kit Composition:

REF		LOT	
KEAA651	LYO CONTROL 1	K354	2027-08
KEAA661	LYO CONTROL 2	K355	2027-08
KEAA511	MTP	K339	2027-09
KEAA611	LYO CAL	K353	2027-08
KEAA751	CONC DILBUF	K338	2027-09
KEAA771	SUBS	K352	2027-08

Valid edition of the instructions for use: 2014-04

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2025-10-27

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl